

# EDK

## Electrostatic Decontamination Kit



## Decon will Never Be the Same...

The Electrostatic Decon Kit (EDK) combines top of the line technology from Victory Innovations and First Line Technology to offer users the easiest, most efficient way to decontaminate threats.

Stored in a military grade hardshell case from SKB Cases, the electrostatic spraying system can easily be charged while staying organized inside the case for rapid deployment when needed.

### MG200 Electrostatic Sprayer

Electrostatic spraying is the process of adding an electric charge to the liquid droplets when they are sprayed. The result is the decon solution will envelop and adhere to the target. As the charged droplets hit the surface they are forced to create an even spread, including hidden and shadowed areas.

### DAHLGREN DECON + FIBERTECT

provides high efficacy decontamination products that can decontaminate chemical warfare agents, toxic industrial chemicals (TICs), toxic industrial materials (TIMs), as well as biological and radiological agents from personnel, personal protective equipment (PPE), sensitive equipment, and critical infrastructure.

### Kit Includes:

- + MG200 Electrostatic Sprayer
- + (x2) Dahlgren Decon
- + (x8) Fibertect Poly Wipes / Pads
- + 1000ml Reservoir
- + 8 hour extended battery
- + 110v Battery Charger
- + Protective Case

## Dahlgren Decon neutralizes 100% of Fentanyl and Carfentanil in under five minutes.

*Testing done on fentanyl Freebase and carfentanil Oxalate at MRI Global, June 2017.*



800.892.1099 | federalresources.com



# EDK Electrostatic Decontamination Kit



## Professional Cordless Electrostatic Sprayer


### Key Highlights:

- + Lightweight and Portable - 3.8lbs.
- + Quiet — Suitable for use in any environment.
- + Removable tank holds 38oz. of liquid.
- + Charge on/off for wetter applications.
- + Save Time, Save Money
- + Increased Performance

### The Electrostatic Difference

Conventional Sprayer

Victory Sprayer



### THE STATS:

pH Level ..... **6.64 (noncorrosive)**  
 Shelf Life ..... **10 Years**  
 Unknown Toxicity ..... **0%**  
 Dwell Time ..... **Under 5 Minutes**  
 Storage ..... **-40° - 160° F**

- + **Non-Corrosive**
- + **Universal**
- + **Neutralization**
- + **Rapid Response**

Dahlgren Decon's pH is truly neutral at 6.64.  
 Dahlgren Decon can be used in conjunction with VICTORY 200  
 Dahlgren Decon neutralizes fentanyl with no residual byproduct  
 In a 1:1 ratio, Dahlgren Decon yields 100% efficacy in less than 2 minutes against fentanyl

	People	Weapons	Sensitive Equipment	Individual Equipment
<b>DAHLGREN DECON</b>		✓	✓	✓
<b>FIBERTECT</b>	✓	✓	✓	✓



800.892.1099 | federalresources.com

